

# Fecal Occult Blood Test Devices

**Pro  
ADVANTAGE**<sup>®</sup>  
by **NDC**

Rapid, convenient, and qualitative methods for asymptomatic occult bleeding. Available in guaiac slide – with both color enhanced sensitivity and readability.

- On board positive and negative monitors
- Includes enhanced color sensitivity and readability
- CLIA waived and FDA 510(k) cleared to market



Item #	Description
<b>P080008</b>	<b>Fecal Occult Blood Patient Pack – Guaiac Slides</b> <i>Includes: 100 Triple Slide Tests, 6x10mL Developers, Applicators, Foil-Lined U.S. Postal Service approved envelopes. 2x50 per box</i>
<b>P080018</b>	<b>Fecal Occult Blood Lab Pack – Guaiac Slides</b> <i>Includes: 100 Single Slides Tests, 2x10mL Developers, Applicators. 100 per box</i>
<b>P080028</b>	<b>Fecal Occult Developer Solution – for Guaiac Slides</b> <i>Includes: 8x10mL Developers. 8 per box</i>

**Distributed By:**



CPT Code: 82270