

& Urine / Serum hCG Pregnancy Test Devices

Pro
ADVANTAGE[®]
by **NDC**

Rapid, one step tests for the qualitative detection of human chorionic gonadotropin (hCG) in serum and or urine. Three variations of Pro Advantage hCG tests are available: Urine Cassette Devices, Urine/Serum Pregnancy Devices and Urine Test Strips.



- >99% Sensitivity; >99% Specificity
- 25 mIU/mL minimum hCG detection level
- Accurate results in 3 minutes for urine and 5 minutes for serum
- CLIA waived for urine and FDA 510(k) cleared to market.



For Urine Samples,
Serum: Moderately Complex
CPT Code:
81025 Urine, 84703 Serum

Item #	Description
--------	-------------

PREGNANCY TEST DEVICES

- | | |
|----------------|--|
| P080063 | Urine hCG Pregnancy Cassette Device
<i>Includes: 25 individually packaged devices and droppers. 25 tests per box</i> |
| P080021 | Urine/Serum hCG Pregnancy Cassette Device
<i>Includes: 25 individually packaged devices and droppers. 25 tests per box</i> |
| P080023 | Urine hCG Pregnancy Test Strips
<i>Includes: 20 Test Strips per canister, 5 canisters per box</i> |

Distributed By: