

A^{Pro}
ADVANTAGE[®]

by **NDC**

Ultrasound and Lubricating

Diagnostic Aids



Ultrasound Transmission Gel

- Thick, no-slip formula eliminates air between the transducer and skin offering exceptional coupling efficiency for more effective treatment and diagnosis.
- Formulated to be acoustically correct for the broad range of frequencies used in medical ultrasound applications.
- Non-irritating hypoallergenic formula will not sting or chafe patient tissue.
- Salt-free and formaldehyde-free.
- Water soluble formula prevents damage to equipment or clothing.

Lubricating Jelly

- Lubricates body orifices to facilitate the entry of diagnostic and therapeutic devices.
- Greaseless, water soluble, and bacteriostatic.
- Available in single use packets and a flip top tube.
- Contents remain sterile in unused, unopened packages.

Item #	Description
	Ultrasound Transmission Gel
P301008	8 oz. Bottle with Flip Cap
P301105	1.3 Gallon (5 L) Collapsible Container with Empty Bottle
	Lubricating Jelly
P903100	3 gm. Packets, Sterile, 144/box
P903200	4 oz. Flip-top Tube, Sterile